41. Betazerfall - Beta Decay

Überblick

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Overview

In this experiment we are considering the β -decay, and different measurements are performed using a β -emitter. In a β -decay, either a proton is converted to a neutron (β^+ -decay), or a neutron is converted to a proton (β^- -decay). The emitted electron (or in the case of the β^+ -decay a positron is emitted) is counted by means of a detector.

The goal of this experiment is threefold:

- By measuring the activity with the detector (events per second), and knowing the geometry of the setup, it is possible to calculate the instantaneous activity of the beta emitter. By knowing the half-life of the decaying material, it is possible to obtain the activity at any given time.
- The absorption curve should be measured. From that curve we can calculate the range of the electrons in the absorbing material (here: aluminium).
- By analysing the absorption curve, the maximal initial energy of the electrons can be calculated in various ways.

1 Basics

1.1 The Structure of the Atomic Nucleus

Over the last century, numerous experiments and theoretical considerations have led to very clear ideas about the structure of atoms. They consist of the Z -fold positively charged *atomic nucleus*, which carries most of the atom's mass, and the Z -fold negatively charged electron shell. The nuclear radius is about 10^5 times smaller than the atomic radius, which is 10^{-14} to 10^{-15} m.

The nucleus is a composite object with an internal structure. It is composed of protons and neutrons (the so-called *nucleons*). The positively charged *proton* and the electrically neutral *neutron* have almost the same rest mass (about 1837 electron masses). The number Z of protons of a nucleus indicates its charge and thus the order number of the relevant chemical element, which is the same as the number of electrons (since the atom itself is electrically neutral). The number of nucleons $A = Z + N$ (with N = number of neutrons) determines the atomic weight. The contribution due to the mass of the electrons is

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very small.

It is known that the nuclear forces can not be of classical nature: The gravitational forces, for example, are far too weak to play a role in the nuclear structure, in spite of the small distances involved. The nuclear force, nowadays referred to as *strong interaction*, has a short range (order of magnitude 10^{-15} m), is attractive for larger distances and strongly repulsive for smaller distances. Accordingly, the nucleon spacing in the nucleus is rarely less than about 0,4 fm (1 fm = 10^{-15} m). The binding energy per nucleon is practically constant, i.e., independent of the position of the corresponding element in the periodic system. It is approximately 8 MeV per nucleon ($1 \text{ MeV} = 10^6 \text{ eV} = 10^6 \cdot 1.602 \times 10^{-19} \text{ J}$).

Apart from the problem of an exact description of the nuclear forces, also the tools for a complete treatment of the *multi-body problem*, such as the treatment of the atomic nucleus, are generally lacking. For this reason it is necessary to rely on *models* (as it is the case in many fields of physics), which lead to correct results for surprisingly many - but not all - aspects. In particular, the so-called *drop model* and the *shell model* of the atomic nucleus are very popular. In the drop model, the stability and dynamics of a complex atomic nucleus is described analogously to the behaviour of a drop of a liquid. The model has proven to be very useful especially in the description of nuclear fission processes. In the shell model, the motion of the individual nucleons in the nucleus is treated analogously to the motion of the individual electrons in an atom. For the atom, as we know, this treatment (together with the Bohr postulates or wave mechanics) leads to different but stable electron configurations, the so-called electron shells or orbitals. In analogy to this, there are also several configurations for the atomic nuclei which are characterized by particular stability. These are nuclei in which either the proton number Z or the neutron number N takes one of the following "magic values":

$$Z \text{ or } N = 2, 6, 8, 14, 20, 28, 50, 82 \text{ and } 126$$

Although it was possible to develop a formalism which describes these "magic values" and many other properties of the atomic nuclei, one can not talk about a quantitatively exact description of the atomic nuclei in the shell model either. It is not within the scope of this manual to treat the nuclear structure in detail. To do so, please refer to standard textbooks. Here, it should just be mentioned that protons and neutrons are no "elementary particles" because they are composed of so-called *quarks*. The original quark model has developed into *quantum-chromodynamics*. This model appears to be capable of explaining the existence and many properties of various particles which were formerly considered being elementary particles (elementary means not composed of other particles).

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The natural and the artificial nuclear instabilities are of great importance when studying the nuclear structure. These include the α -, β -, and γ -radiation. In the natural α -radioactivity, which is observed for $Z > 82$ (Pb), an α -particle (which is the nucleus of a helium atom, i.e., two protons and two neutrons) is emitted. The doubly charged α -particles have a very short range (few cm in air), but a large effect when hitting a surface. In the β -decay, a neutron (proton) is converted to a proton (neutron), under the emission of an electron (positron), depending on the type (β^- or β^+). The β^- -decay will be examined more closely in this experiment. The γ -radiation describes the emission of electromagnetic waves usually due to de-excitation of (decay-) nuclei, with energies above some 100 keVs.

Another important method for studying nuclear structures is provided by the *nuclear reactions*. By bombarding nuclei with neutrons, protons, and other particles, transformations can be induced and observed. If the observations are carried out systematically, information about the reaction energy of a process and the quantized energy states in the nuclei can be inferred.

Probably the most important and also the best investigated type of nuclear instability, however, is the β -decay. It is described in the framework of the so-called *weak interaction*. In the β -decay, the unstable nuclei emit a particle, which can be identified as fast electron. In other cases, a positive electron, the so-called *positron*, is observed. The positron has the same charge e (though with the opposite sign!) and the same absolute value of the e/m ratio and spin as the electron. The positron is the so-called "anti-particle" of the electron. Its existence was predicted by P. Dirac (1928) on the basis of a relativistic extension of wave mechanics.


Fig. 1: Diskretes Energiespektrum des α -Zerfalls. / Discrete energy spectrum of the α -decay.


Fig. 2: Kontinuierliches Energiespektrum des β -Zerfalls. / Continuous energy spectrum of the β -decay.

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The energy of the atomic nuclei (like that of the electrons in the atoms, see also experiment “24 – Spectroscopy”) is not continuous, but it only assumes discrete values (quantization). Thus, the energy spectrum of the emitted α -particles in a radioactive decay process shows the character of a *line spectrum* (Fig. 1).

When first discovered, it was astonishing that the β -decay particles did not show a uniform (discrete) energy spectrum. The observed spectrum was continuous with a sharp upper *energy limit* E_{max} , as shown in Fig. 2. This was first discovered by J. Chadwick in 1914. Such a behaviour was all the less understandable, as it could also be demonstrated experimentally that an identical, well-defined energy quantity is transferred in each individual decay process. It actually seemed as if the laws of the conservation of energy (and momentum) were violated during the β -decay. W. Pauli postulated (1932) that a second, uncharged particle has to be emitted during the β -decay. This particle was later referred to as a neutrino, denoted ν .

In 1934, E. Fermi calculated the expected energy spectrum of the β -decay from Pauli’s neutrino hypothesis, and obtained good agreement with the experiment. This could be interpreted as an indirect proof of the existence of the electron neutrino ν_e (there exists also the muon neutrino, ν_μ , and a tau neutrino, ν_τ). However, the direct experimental detection of the neutrino could not be achieved until 1955 by F. Reines and C. Cowan, because neutrinos are very weakly interacting and thus very difficult to detect¹.

The important β -decay with all involved partners is summarized below. Note that a bar on top of a particle’s name, such as $\bar{\nu}$, denotes the antiparticle. Only for the antiparticle of the electron exists a dedicated name and symbol, i.e., the positron e^+ .

$$\begin{aligned}\beta^- \text{-decay} : n &\rightarrow p + e^- + \bar{\nu}_e \\ \text{on quark-level} : d &\rightarrow u + e^- + \bar{\nu}_e \\ \beta^+ \text{-decay} : p &\rightarrow n + e^+ + \nu_e \\ \text{on quark-level} : u &\rightarrow d + e^+ + \nu_e\end{aligned}$$

¹Dedicated neutrino experiments, such as for example the Super-Kamiokande experiment in Japan, require huge detectors consisting of 50'000 tons of ultra pure water and more than 12'000 photomultiplier. These vast numbers give a hint on how difficult it is to precisely measure properties of neutrinos.

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1.3 The Decay Mode of ^{90}Sr

In this experiment, we investigate the β^- -decay of ^{90}Sr . The isotope ^{90}Sr decays via β^- -decay into ^{90}Y . The maximum decay energy is 546 keV, and the half-life is about 28.8 years. However, ^{90}Y is not stable, but decays further into ^{90}Zr . Zirconium is stable and there are no further decays. The half-life of ^{90}Y is 64 hours, and the branching ratio for the β^- -decay is almost 100%. Only 0.01% of the decays from ^{90}Y are gamma-decays. Strontium and Yttrium are in equilibrium within the source, and per one decay of Strontium we observe one decay of Yttrium as well. The source is covered with a thin foil made from steel, in order to prevent the radioactive material coming out of the source.

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The energy spectrum of the source used in this experiment is given in Fig. 3 (image modified from [3]). The β -particles (i.e., electrons) from the decay of ^{90}Sr have a continuous energy spectrum, and a fraction of these electrons have not sufficient energy to penetrate the thin layer of steel covering our source. The energy of the β -particles from the decay of ^{90}Y

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have much higher energies (although also continuous), and most of them manage to leave the source. Thus, we primarily observe in this experiment the decay electrons of $^{90}\text{Y} \rightarrow ^{90}\text{Zr}$. The maximal decay energy of the electrons from the decay of ^{90}Y , E_{max} , should be determined in task 5.

1.4 Range of the Particles

The energetic electrons originating from radioactive decay processes lose their energy gradually as a result of a large number of collisions with other atoms. The distance they travel until they come to rest is called the *range* R of the electrons.

The range depends on the material itself (e.g. air or aluminium) and the initial energy of the electrons. The range can therefore be used to determine the initial energy of the electrons. In gases, R is inversely proportional to the density ρ . Thus, the product

$$\rho R = f(E) \quad (1)$$

is a function of the initial energy.

Although in solids the density ρ is constant, it has become natural to specify the product ρR , and not only the range R . If the function

$$E = f^{-1}(\rho R) \quad (2)$$

of the absorber material and its density are known, then the initial energy of the electrons can be determined from the range in the material. Fig. 4 shows the function (2) for aluminium.

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In our case, we are mainly interested in the *limit energy* E_{\max} (Fig. 2) of the electrons emitted during the β -decay. It is therefore necessary to experimentally determine the *maximum* range of these electrons.

We denote by N_0 the number of electrons that are emitted into an absorber material and by $N(x)$ the number which pass through a layer of thickness x . For a β -spectrum (with an energy distribution as shown in Fig. 2) it applies:

$$N(x) = N_0 e^{-\mu x} \quad (3)$$

The coefficient μ depends on the maximum energy E_{\max} of the β -spectrum and on the density of the absorber material. They are related by the following empirical relationship:

$$\mu = 17 \frac{\rho}{(E_{\max})^{1.43}} \quad (4)$$

with μ in cm^{-1} , E_{\max} in MeV, and ρ in g cm^{-3} .

The experimental determination of μ provides a way to calculate E_{\max} . If the absorber thickness x approaches the maximum range R_{\max} , a deviation from the exponential behaviour (Eq. 3) occurs, as indicated in Fig. 5. Here, the absorption curve can be approximated by the ansatz

$$N(x) = \text{const.} [E_{\max} - E(x)]^n. \quad (5)$$

$E(x)$ corresponds to the expression (2) (with $R = x$) and, for aluminium, can be extracted from Fig. 4. The exponent n depends on the investigated β -spectrum and in general attains values between 3 and 4. Using equation 5 allows to extrapolate E_{\max} from an experimentally measured absorption curve.

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1.5 Detectors

The detection of radioactive particles is based on the fact that they lose their energy in ionization and

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excitation processes when passing through matter. Common measuring instruments, which can detect individual particles, are the *counter tube* and the *scintillation counter*. Other important detector types are the photographic plate, the Cerenkov counter, the (Wilson) cloud chamber, the spark chamber, semiconductor detectors and superconducting detectors. However, in the following, we will restrict the discussion to the counter tube, as it is used in this experiment.

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A standard counter tube (by H. Geiger 1913, H. Geiger and W. Müller 1928) consists of a cylindrical capacitor (Fig. 6). A straight metallic cylinder serves as a cathode and an electrically insulated, axially extending wire serves as an anode. The housing is hermetically sealed and contains a filling gas. A typical mixture is 90 Torr Argon and 10 Torr Ethyl Alcohol. One end of the cylinder is sealed by a thin entrance window that consists of low-Z (little scattering) material such as Mylar.

If an energetic particle then passes through the window into the interior of the counting tube, it will ionize the molecules of the filling gas. Depending on the magnitude of the voltage V applied to the counter tube, different reactions of the filling gas are occurring. It can be distinguished between four voltage ranges:

1. At very low voltages, most ions *recombine* before they can reach the electrodes. The counter tube practically does not respond at all.
2. At slightly higher voltages, most of the originally generated ions can reach the electrodes. The charge transported thereby provides a measure for the number of incident particles. The counter tube thus acts as an ionization chamber. However, the ion current emitted by a single particle is too small to be detected.
3. If the voltage is increased further, the electric field strength in the vicinity of the wire can become so large that the electrons gain sufficient energy during one free path length, such that they ionize further atoms of the filling gas. Each

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initially created ion can then lead to a localized ionic avalanche („Townsend discharge“). There is so-called *gas amplification*. The discharge breaks off after about 10^{-6} seconds. As long as the gas amplification factor A (number of secondary ions generated per primary ion) is independent of the overall number of primary ions, the transported charge is proportional to the number of initially created ions. This is called the *proportionality range*.

4. By increasing the applied voltage further, we enter the so-called *Geiger plateau* range of the counter tube. Each discharge is then no longer localized, but propagates over the entire counting tube. Due to the special composition of the filling gas and the geometry of the counting tube, it is possible to suppress autonomous discharges (i.e., no discharge if no particle enters the counter tube). Each discharge stops after about 0.1 millisecond.

At this point, there is always approximately the same amount of charge flowing through the counter tube, irrespective of the number of primarily produced ions. In order to detect the signal, the short voltage drop over the resistor R (see Fig. 6) is measured rather effortlessly with an *electronic counter*. Since the counting tube can only register a new particle after a time of approximately 10^{-4} seconds, the maximum count rate obtainable in the Geiger plateau is also of the order of 10^4 s^{-1} .

The Geiger plateau begins at voltages of the order of 1000 V and is a few 100 V wide. It is expressed in the so-called *counter tube characteristic* (pulse rate with constant irradiation as a function of the voltage) through a correspondingly wide plateau. The small slope of some % per 100 V is due to the fact that as the voltage increases, an ever greater part of the counter tube volume becomes sensitive. If the voltage is increased even further, self-discharging occurs.

1.6 Different Units of Radiation and Biological Effects

Radiation which impinges on biological tissue is in general problematic and usually harms the tissue. The reason for this is that radiation loses energy in the tissue, and the deposited energy can affect the cell division or damage the DNA, such that e.g. cancer can develop. In order to quantify the effects of radiation on biological tissue, it is not sufficient to only indicate the number of incident particles, i.e., the activity, but other units are used, which is discussed in what follows.

The activity, which is the actual number of decays in a certain time interval, is measured in „Becquerel“ (Bq), which is number of decays per second:

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Another common unit is the outdated unit “Curie” (Ci), with $1 \text{ Ci} = 3.7 \cdot 10^{10} \text{ Bq}$. A simple counting detector, such as the one used in the present experiment, cannot distinguish between the different types of radiation. As long as the radiation fulfils the detector’s requirement to be detected (i.e., a certain minimal energy, a certain energy deposition within the detector), the detector counts it.

But a unit which is only counting radioactive decays does not make a statement about the energy which has been deposited in a certain volume. Thus, when quantifying the deposited energy, the unit “Gray” (Gy) is used, with $1 \text{ Gy} = 1 \text{ J/kg}$. But because different radiation has different means of depositing energy, and thus different effects on (biological) material, the unit which actually matters when it comes to biological effects, is “Sievert” (Si). The Sievert is composed from the energy dose, weighted by a number quantifying the effect of the corresponding particle, which is called the “radiation weighting factor”, usually denoted W_R .

The radiation weighting factor for x-rays, gamma rays, beta particles and muons is 1. For neutrons, it strongly depends on their energy, and varies between 5 and 20. For protons we have $W_R = 5$, and for alpha particles and other heavy nuclei $W_R = 20$.

To simplify the calculations, the so-called “dose factor” takes all these factors into account, and additionally also the way of obtaining the radiation (e.g. oral incorporation, inhalation etc.). Additionally, because each nuclide emits radiation at different energies, the dose factor provides the conversion factor from Bq to Sv for each nuclide². Some example dose factors are listed in table 2. Note the huge difference between the β -emitters and americium, which is an α -emitter. Even without working with artificial radioactive sources, we are exposed to radiation everyday. The most dominant source of radiation is radon (ca. 3.2 mSv/year), followed by radiation due to medical reasons (1.2 mSv/year). In total, the average annual dose in Switzerland is 5.5 mSv per person [2].

²There are even more detailed factors, taking into account e.g. specific tissues or organs. However, for practical purposes, it is usually sufficient to look up the dose factors in corresponding tables.

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2 Tasks

In this experiment, a commercial ^{90}Sr -source is used. It is encapsulated in a small metal cylinder, covered with 0.1 mm stainless steel. The source is then mounted in an acrylic glass shielding, below the counter tube (see also Fig. 7). The activity of the source should be determined in one of the tasks, but as reference value it is between 50 and 100 kBq.

1. Measure the counter tube characteristic and determine the working point voltage (Geiger plateau) V_S as well as the plateau slope $\Delta N/\Delta V$.
2. Determine the current activity of the radioactive source.
3. Calculate the time it takes for the activity to drop to 10 kBq, calculate what the activity was 70 years ago, and how long you have to be fully exposed to the pure, naked source to reach the average annual dose in Switzerland (assuming incorporation).
4. Measure the count rate N as a function of the absorber thickness x and determine the absorption curve.
5. From the absorption curve, determine the maximum β -energy of ^{90}Y in different ways.

Attention: The sample used here is a *closed* source with a thin layer of metal covering the radioactive material. Nevertheless, the sample must never be touched or removed from the shielding. In the event of (possible) contamination, inform the assistant or the helpdesk *immediately*, *without* further spreading the contamination!

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3 Conduction of the Experiment

The measurements are carried out in the closed container. The setup is shown in Fig. 7. The container must *not* be opened during the experiment! Absorber sheets made of aluminium of different thickness can be placed on a slider. The absorber is positioned correctly if one end of the slider coincides with the surface of the container.

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3.1 Counter Tube Characteristic

The slider, which later on serves to hold the absorbers, acts as an aperture. The number of counts $N(V)$, i.e. the number of registered radioactive decays is measured as a function of the applied voltage (circuit according to Fig. 6). The measurement time t should be chosen in order to obtain reasonable statistics (e.g. 1% accuracy, see also section 4). Start with $V = 200$ V and increase the voltage in steps of 100 V. The voltage is measured with an electrostatic voltmeter. Plot the number of detected counts as a function of the applied voltage. For our measurements, the counter tube must be operated in the Geiger plateau. The working point voltage V_{ideal} should therefore be chosen to be at the centre of the plateau.

Plot the count rate N as a function of the applied voltage V and fit a straight line in the plateau region. Also, indicate the voltage dependence of the count rate $\Delta N/\Delta V$.

Does a small $\Delta N/\Delta V$ yield a better or a worse detector?

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$$\Omega = 4\pi \sin^2\left(\frac{\theta}{2}\right).$$

3.2 Activity of the Source

The term *activity* of a radioactive source describes the number of decays per second. The units of the activity is *Becquerel* (Bq). A source has an activity of 1 Bq, if there is 1 decay per second.

In order to determine the activity of the sample, we first have to determine the effective activity N_{eff} of the β -particles which enter the detector. To do so, the total counts in a certain time interval, e.g. 60 seconds, is recorded, and we call it N . Then, the thick 8 mm steel slider is used to shield the detector completely from the β -radiation. The detector now only counts background radiation. We denote the background count rate with N_{BG} . We further assume that every particle hitting the detector is registered, thus the responsivity of the detector is 100%. The effective activity N_{eff} of the source is then $N_{\text{eff}} = N - N_{\text{BG}}$. We assume that the radioactive source radiates evenly in all directions in space. However, the detector only detects a fraction of the radiation because it cannot cover all directions around the source. We call this fraction of the covered directions the acceptance ε of the detector. The acceptance is between zero and one, $0 \leq \varepsilon \leq 1$.

To calculate the acceptance we can use the concept of *solid angle* Ω : The solid angle Ω is the subarea A of a sphere surface with radius r , divided by r^2 :

$$\Omega = \frac{A}{r^2}$$

The full solid angle Ω_{max} (i.e. all directions) is $\Omega_{\text{max}} = 4\pi$. The solid angle Ω of a cone with half angle θ (see figure 8) is given by

$$\Omega = 4\pi \sin^2\left(\frac{\theta}{2}\right).$$

Since the source radiates evenly in all spatial directions, the acceptance ε is given by $\varepsilon = \Omega/4\pi = \sin^2(\theta/2)$.

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The point-like source is $r = 101$ mm away from the aperture of the steel slider (see also figure 7. Measure the radius r of the aperture in the slider in order to determine the half angle of the cone θ and the acceptance ε .

Together with N_{eff} , this value allows you to calculate the current activity of the source. Note, however, that some electrons do not penetrate the thin steel layer in front of the source. We correct for this as described below:

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In Fig. 9, the electron energy loss in 0.1 mm stainless steel is given as function of electron energy. As can be seen, it is a good approximation that all electrons with energies below 0.2 MeV cannot penetrate the thin steel window. From simulations³, this amounts to about 55% of all electrons emitted by the decay of ^{90}Sr that do not reach the detector, which - from geometrical considerations - would have reached the detector if the foil was not present. Analogously, 10% of all electrons emitted by the decay of ^{90}Y do not reach the detector, although their geometric trajectories would meet the detector.

Since the half-life of ^{90}Sr is many times longer than that of ^{90}Y , we can speak of a double decay. If a strontium atom decays, an yttrium atom decays basically immediately afterwards. So if the source emits 100 beta particles, an average of 50 would come from a strontium decay and 50 from an yttrium decay.

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³From RADAR-data [4], the energy spectrum of ^{90}Sr and ^{90}Y was calculated. Then, the fraction of electrons with energies below 0.2 MeV was approximated by integration.

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3.3 Activity and Dose Calculations

Radioactive decay processes are exponential in character. The number of atoms $n(t)$ decreases exponentially $n(t) = n_0 e^{-t/\tau}$. We call τ the lifetime, or $1/e$ -time, because after a time $t = \tau$ there are still n_0/e atoms present.

- a) Calculate the $1/e$ times of Sr and Y based on the half-lives.

In exponential processes the change of the number of atoms $\partial n/\partial t$, and thus the activity, is proportional to the number of atoms, $\partial n/\partial t \propto n$. The activity $A(t)$ therefore also decreases exponentially $A(t) = A_0 e^{-t/\tau}$.

- b) Based on the measured current activity A_0 of the radioactive source, calculate how long it will take until the activity has decreased to $A = 10$ kBq.

Even when we do not come into direct contact with radioactive sources, we are constantly exposed to a certain amount of background radiation. This radiation comes from cosmic sources and radioactive elements in the air, such as radon. In Switzerland, the average annual radiation dose is $S = 5.5$ mSv.

- c) Calculate how long it takes to reach the annual radiation dose of 5.5 mSv if the radioactive source of the experiment is swallowed (i.e., incorporated) completely without the steel film.

Use the dose factors from table 2. It can be assumed that the dose coefficients have the unit Sv/count instead of Sv/Bq . Since ^{90}Y has a much shorter half-life than ^{90}Sr , it can be assumed that the activities are identical.

3.4 Absorption Curve

The maximum energy of the beta particles, E_{\max} , shall be determined in this task. The higher the energy of the particles, the thicker layers they can penetrate. Figure 4 shows how much energy is needed to penetrate a certain material thickness. More precisely: Since the density of the material is also of great importance for the range of the particles in this material, the range is given as a function of the distance-density product.

Place aluminum plates of different thicknesses x in the steel slider and measure the count rate N over a period of time, for example $t = 60$ s. Then calculate for each measurement the product density times thickness ρx , where the density of aluminum $\rho = 2.69 \text{ g cm}^{-3}$. The thickness of the aluminum platelets is measured with the caliper. Subtract the background from the measured count rate N to obtain N_{eff} and calculate the effective activity. From

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figure 4 you should also read the energy loss $E(x)$ that the electrons experience when passing through the respective absorber platelets.

After you have made measurements for different platelet thicknesses x , you should display the results graphically. To do this, make a plot with the logarithm of the activity $\ln(N/t)$ on the y-axis and the plate thickness x on the x-axis. We expect the activity to decrease at first, because electrons with less energy do not reach the detector, until finally the activity flattens and equals the background value, as shown schematically in figure 10.

- a) Determine from eye at which thickness x_{max} the curve flattens. The energy at this thickness corresponds to the maximum energy of the electrons E_{max} .

When penetrating the steel foil around the radioactive source, the electrons already lose energy.

- b) Determine this energy loss using Figure 9.
- c) Calculate the maximum energy of the electrons directly in source $E_{\text{max source}}$.

Now, fit the first few values (where the activity was above the background, i.e. those with $x < x_{\text{max}}$) with a straight line. The slope of this line corresponds to μ from equation 3, $N(x) = N_0 e^{-\mu x}$.

- d) Enter the slope μ and calculate the maximum energy of the electrons E_{max} using equation 4, $\mu = 17\rho/(E_{\text{max}})^{1.43}$. Again, we measure the energy after the steel foil surrounding the source.
- e) Calculate the maximum energy at the source based on the measurement of the slope.

We have now determined the maximum energy of the electrons by two independent methods.

3.5 Determining E_{max} of ^{90}Y for physics students

A better result can be achieved by using equation 5. First, we subtract the background: Subtract N_{BG} from N : $N_{\text{eff}} = N - N_{BG}$. Plot N_{eff} as a function of $E_{\text{max}} - E$ in a double logarithmic scale (see sketch in Fig. 11 (left)). For E_{max} we use the value obtained before (method 1).

From the average slope, determine n in equation 5. Plotting $(N_{\text{eff}})^{1/n}$ as a function of E (Fig. 11 (right)) yields a better value for E_{max} (which is the intersection of the interpolation line with the abscissa).

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4 Statistical considerations

The measurements of the count rate (number of decays per second) will become more and more accurate the longer the measurement lasts, i.e., the larger the number of measured counts is. What is now the uncertainty for such a measurement?

Assume that $\langle n \rangle$ is the count rate, averaged over a long time. If we count the detected events several times during t seconds, then the measured count rate will fluctuate around the true value. This fluctuation is described by the *Poisson*-distribution. A dedicated discussion of the Poisson distribution can be found e.g. in [5].

For the Poisson-distribution, the average uncertainty $\langle \Delta n \rangle$ (mean value of many identical measurements) is:

$$\langle \Delta n \rangle = \sqrt{n}$$

where n is the number of counts which are detected during the time interval t .

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The relative uncertainty will be

$$\frac{\langle \Delta n \rangle}{n} = \frac{1}{\sqrt{n}}.$$

Thus, in order to reach a precision of 1 %, we should count up to $n \approx 10^4$ events.

5 Questions for physics students

1. Describe the various ranges of operation of a gas detector. Which voltage range should be chosen to operate the detector as ionisation chamber, proportional chamber or Geiger-Mueller counter tube?
2. What determines the energy resolution of the detector? A particularly good energy resolution can be achieved with semiconductor detectors. Why? How do scintillation detector work?
3. What other important aspects of detectors can you think of?
4. How can spatial resolution be achieved with counting detectors?
5. Beta particles are either negatively charged (=electrons) or positively charged (=positrons). How are positrons used nowadays to investigate and examine materials?
6. Are neutrons also subject to β -decays?
7. What is the difference between γ -rays and X-rays?

Literatur

- [1] SCHWEIZERISCHE STRAHLENSCHUTZVERORDNUNG DES BUNDES (StSV), *Version vom 26. April 2017, Stand 1. Februar 2019*
online: <https://www.admin.ch/opc/de/official-compilation/2017/4261.pdf>; Zugriff am 10. Januar 2020
- [2] HINTERGRUNDINFORMATION: ORTSDOSISLEISTUNG UND DOSIS, *Version vom 21. August 2015*
online: https://www.naz.ch/de/downloads/Factsheet_ODL_Dosis_de.pdf; Download am 10. Januar 2020
- [3] A SYSTEM FOR THE MEASUREMENT OF ELECTRON STOPPING POWERS: PROOF OF PRINCIPLE USING A PURE β -EMITTING SOURCE, *Radiation Physics and Chemistry*, **149**, 134-141 (2018), TIMOTHY ROY, FREDERIC TESSIER, MALCOM MCEWEN
- [4] BetaSpec.zip. DOWNLOAD AM 16. JANUAR 2020, von RADIATION DOSE ASSESSMENT RESOURCE (RADAR); <http://www.doseinfo-radar.com/RADARDecay.html>
- [5] PRACTICAL PHYSICS, G. L. Squires, 1985, (3rd edn). Cambridge: Cambridge University press.